

Specimen Collection & Transport Kits for Drugs of Abuse Testing

- Features Wide-Mouth Collection Cup
- Includes 45 ml Tamper Evident Vial
- Single and Split Kit Configurations
- Bag and Box Kits Available
- Customization available upon request

Therapak assembles Non-DOT, DOT and SAMHSA compliant collection and transport kits using the popular 45 ml capacity tamper evident vial. The vial is constructed of virgin polypropylene that results in excellent clarity for visibility into the collected specimen container. Each vial is graduated at the 15 ml and 30 ml capacity level. Kit configurations include a two-pouch specimen bag with absorbent, a 180 ml beaker style collection cup with pour spout and either one or two 45 ml vials depending on whether your requirements call for single or split collection kits. Kit components are heat-sealed in a bag to protect the integrity of the components.

Individual single and split collection box kits are also available. To customize the box or specimen bag with logo information, contact us at info@therapak.com to review minimum volume and lead times. All kits are manufactured in FDA-registered facilities in either California or Georgia in accordance with cGMP-compliant processes.



#39467 Single 45 ml Vial Bag Kit



#39468 Split 45 ml Vial Bag Kit

Ordering Information

| Part No. | Description | Packed |
|----------|---------------------------|----------|
| 39467 | Single 45 ml Vial Bag Kit | 100/case |
| 39468 | Split 45 ml Vial Bag Kit | 100/case |
| 39469 | Single 45 ml Vial Box Kit | 48/case |
| 39470 | Split 45 ml Vial Box Kit | 48/case |

Contact Therapak to learn more about online ordering, requisition printing and direct ship programs.

info@therapak.com • (888) 505-7377 • www.therapak.com

Bulk Supply Kits including 45ml Tamper Evident Vial

- Reduces outbound freight costs to client locations
- Nested collection cups individually wrapped with temperature strip
- Order Multi-sample shipping boxes separately for shipping specimens back to the laboratory

Therapak's Bulk Supply Kits provide the most economical drug testing supply solution in the industry. Bulk packaging means that items are nested and packaged in a manner to reduce assembly costs while providing case quantity dimension reductions as compared to other kits being offered in the substance abuse testing market. Freight costs are a major consideration in transporting kits to your facility and from your facility to client locations. These space-saving Bulk Supply Kits are the answer to your cost reduction initiatives. Choose from three different packaging configurations for both single and split collection applications. Packaging configurations include case quantities of 25, 50 and 100.



Single Bulk Supply Kits include individually wrapped beaker collection cups with temperature strips, specimen transport bags with absorbent and 45ml tamper evident vials. Split Bulk Supply Kits include individually wrapped beaker collection cups with temperature strips, specimen transport bags with absorbent and two 45ml tamper evident vials. For more information, contact us at info@therapak.com to learn about customized kits, requisition printing and direct shipping via Therapak's Online ordering systems.

| Part No. | Description | Packed |
|----------|------------------------------|----------|
| 40054 | Single 45 ml Bulk Supply Kit | 25/case |
| 40055 | Single 45 ml Bulk Supply Kit | 50/case |
| 40056 | Single 45 ml Bulk Supply Kit | 100/case |
| 40057 | Split 45 ml Bulk Supply Kit | 25/case |
| 40058 | Split 45 ml Bulk Supply Kit | 50/case |
| 40059 | Split 45 ml Bulk Supply Kit | 100/case |

Single Specimen Collection Kits

- Kits include a 95kPa certified 90 ml container
- Kits are individually bagged
- Kits with transport bags are constructed with two sealable, tamper evident pouches and absorbent material inside the specimen pouch

Therapak drugs of abuse collection and transport kits are available in various configurations. Kits can be used in conjunction with on-site testing products or remote site collection when specimen is transported to a laboratory for analysis.



| Part No. | Description | Packed |
|----------|--|----------|
| 40012 | 90 ml container kit with temperature strip | 100/case |
| 40003 | 90 ml container kit with evidence seal | 100/case |
| 40004 | 90 ml container kit with temperature strip and evidence seal | 100/case |
| 40005 | 90 ml container kit with evidence seal, temperature strip, & transport bag | 100/case |
| 40010 | 90 ml container kit with temperature strip and transport bag | 100/case |

Split Specimen Collection Kits

For drugs of abuse testing applications requiring two samples to be collected from a single donor, Therapak's split specimen collection and transport kits are available in either a split bag kit (40002) or the same kit within a transport box (40011).

Kits include: 6.5oz beaker style collection cup with pour-spout and temperature strip, 2 - 60 ml containers with integrity seal, and a transport bag with absorbent.



| Part No. | Description | Packed |
|----------|-------------------------------|---------|
| 40002 | Split Collection Bag Kit | 50/case |
| 40011 | Split Collection Box Kit, Box | 48/case |

Blood-Alcohol Collection & Transport Kit

Therapak's Blood-Alcohol kit includes all necessary components for the collection and transport of blood samples collected from suspect DUI offenders. Kit contents includes 2 - 16x100mm, 10 ml gray top evacuated blood tubes (Sodium Fluoride - Potassium Oxalate preservative in each tube), blood draw needle, needle holder, iodine prep pad, plastic tube holder, 2-evidence seals, box seal, transport bag with absorbent material, marking labels, and transport box. All components are sealed in a bag inside the box. Compliant for shipping by ground or air in accordance with Exempt Human Specimen shipping regulations.



| Part No. | Description | Packed |
|----------|--|---------|
| 38107 | Blood-Alcohol Collection & Transport Kit | 24/case |

Individual Supplies for Drugs of Abuse Testing

| Part No. | Description | Packed |
|----------|-------------|--------|
|----------|-------------|--------|

| | | |
|-------|---------------------------|----------|
| 74829 | 45 ml Tamper Evident Vial | 700/case |
|-------|---------------------------|----------|



| | | |
|-------|---------------------------|----------|
| 74811 | 90 ml Screw Cap Container | 300/case |
|-------|---------------------------|----------|

| | | |
|-------|---|----------|
| 39415 | 90 ml Screw Cap Container with temp strip | 300/case |
|-------|---|----------|



| | | |
|-------|-----------------------|----------|
| 74850 | 180 ml Collection Cup | 500/case |
|-------|-----------------------|----------|

| | | |
|-------|---------------------------------------|----------|
| 39408 | 180 ml Collection Cup with temp strip | 500/case |
|-------|---------------------------------------|----------|

| | | |
|-------|--|----------|
| 39473 | 180 ml Collection Cup, individually wrapped, with temp strip | 500/case |
|-------|--|----------|



| | | |
|-------|---------------------------------------|------------|
| 15610 | Specimen Transport Bag with absorbent | 2,000/case |
|-------|---------------------------------------|------------|



| | | |
|-------|--|----------|
| 39554 | Tamper Evident Seal, Red 3/4" x 6-5/8" | 100/case |
|-------|--|----------|

| | | |
|-------|--------------------------------------|----------|
| 40008 | Tamper Evident Seal, Red 1" x 6-5/8" | 500/case |
|-------|--------------------------------------|----------|



| | | |
|-------|-------------------------|------------|
| 10300 | Absorbent Sheet, 3 x 4" | 1,000/case |
|-------|-------------------------|------------|



| | | |
|-------|---|----------|
| 40006 | Temperature Monitoring Strip, 90-100° F Includes Celsius and Fahrenheit readings | 100/case |
|-------|---|----------|



| | | |
|-------|--------------------------|-----------|
| 40060 | Multi-Sample Shipper Box | 50/bundle |
|-------|--------------------------|-----------|





Sample Cups for Hitachi Chemistry Analyzers

Free-standing sample cups are precision molded of virgin clear polystyrene to exacting OEM specifications and tolerances for use with Hitachi chemistry analyzers. Conical bottom ensures maximum sample recovery of serum sample by pick-up probe. For use with Hitachi 704, 705, 706, 712, 717, 736, 737, 747, 911, and ES-300 Chemistry Analyzers.

| Part No. | Description | Packed |
|----------|-----------------------------------|----------|
| 31507 | Sample Cups for Hitachi Analyzers | 5,000/cs |

Sample Cups for Roche Chemistry Analyzers

Sample cups are designed for use with RocheCOBAS BIO, FARA, and MIRA chemistry analyzers. Precision molded, of virgin high-density polyethylene, to exacting OEM specifications and tolerances. Cups feature conical wells designed for complete sample pickup, and attached caps to minimize sample evaporation. Integral cap has thin membrane section for easy penetration by sample probe. Five available colors allow for quick sample identification of routine chemistries, stats and controls.



| Part No. | Description | Packed |
|----------|--|------------|
| A-CN130 | Sample Cups for Roche Chemistry Analyzers, natural | 5,000/case |
| A-CP130 | Sample Cups for Roche Chemistry Analyzers, pink | 5,000/case |
| A-CY130 | Sample Cups for Roche Chemistry Analyzers, yellow | 5,000/case |
| A-CG130 | Sample Cups for Roche Chemistry Analyzers, green | 5,000/case |
| A-CB130 | Sample Cups for Roche Chemistry Analyzers, blue | 5,000/case |

Tubes for Toxicology Sample Transfer and Testing

Therapak's range of polypropylene (PP) and polystyrene (PS) sample transfer, storage and testing tubes are available in various sizes (12x75mm to 16x125mm) to enhance drugs of abuse testing workflow. All caps listed are for use with the tube directly above and are made of polyethylene (PE).



| Part No. | Description | Material | Packed |
|----------|------------------------------------|----------|------------|
| 80542 | 12 x 75 mm round bottom test tube | PS | 1,000/case |
| 80547 | 12 mm push cap, natural | PE | 1,000/case |
| 80623 | 13 x 75 mm round bottom test tube | PP | 2,000/case |
| 80557 | 16 x 75 mm round bottom test tube | PP | 1,000/case |
| 94504 | 16 mm push cap, yellow | PE | 1,000/case |
| 80531 | 16 x 92 mm round bottom test tube | PP | 1,000/case |
| 80532 | 16 mm push cap, white | PE | 1,000/case |
| 94626 | 16 x 100 mm round bottom test tube | PS | 1,000/case |
| 80625 | 16 x 125 mm round bottom test tube | PP | 1,000/case |

CONTACT INFORMATION

Los Angeles

1801 Highland Ave, Unit L
Duarte, California 91010

Ph: (626) 357-5900
Fax: (626) 357-5911

Atlanta

4305 Hamilton Mill Road, Suite 200
Buford, Georgia 30518

Ph: (770) 614-2931
Fax: (770) 614-2935

London

Unit 1 Waterway Business Park
Rigby Lane
Hayes, Middlesex UB3 1EY
United Kingdom

Ph: 020 8589 0888
Fax: 020 8589 0850

www.therapak.com • (888) 505-7377