

IDEAL FOR DIAGNOSTIC SPECIMEN SHIPMENTS

Aqui-Pak™ Absorbent Pouch and Tube Separator



Aqui-Pak 4 bay pouch #10314

- Integral part of meeting DOT/IATA Regulations
- Place between primary and secondary container
- Insert up to six (6) tubes in one Aqui-Pak™
- Absorbs up to 150mL of liquid specimen
- Roll it up and place it in a bag or canister.

Aqui-Pak™ was designed for the safe transport of laboratory specimen tubes. Aqui-Pak consists of a multi-layered polymeric material that has extremely high absorption properties. It is very convenient when used in combination with a specimen bag or canister to hold up to six tubes.

Aqui-Pak is an integral solution for packaging Diagnostic Specimens according to the recently published DOT Hazardous Materials Regulations (HMRs) and IATA Dangerous Goods Regulations. Both the DOT and IATA have established packaging requirements for Diagnostic Specimens. Among other instructions for preparing blood and urine tubes (primary containers) inside of transport bags, canisters (secondary watertight receptacles) and boxes (rigid outer container), blood and urine tubes considered Diagnostic Specimens must be packaged with some form of absorbent material that not only separates multiple primary containers to prevent breakage, but cushions them within the packaging.

According to the regulations, absorbent material must be placed between the primary and secondary container and must be capable of absorbing the fluid content of all primary containers within the secondary packaging. Aqui-Pak absorbs, separates and cushions which eliminates the need for multiple product selection for packaging and transporting Diagnostic Specimens. Aqui-Pak is the cost efficient solution for compliance.



Aqui-Pak 6 bay pouch #10316

Part No.	Item	Packaging
10310	Aqui-Pak™ Single bay absorbent pouch	3,000/case
10335	Aqui-Pak™ 2 bay absorbent pouch	4,000/case
10314	Aqui-Pak™ 4 bay absorbent pouch	900/case
10316	Aqui-Pak™ 6 bay absorbent pouch	900/case

FOR MORE INFORMATION

(909) 267-2000/(770) 614-2931

info@therapak.com

Aqui-Pak™ is a trademark of Therapak, LLC.

www.therapak.com