

Dual Click-Tite™ Containers

Therapak's Dual Click-Tite™ Containers are ideal for maintaining the integrity of biological specimens during transit and storage.

Durable Design

Composed of 100% virgin FDA approved high clarity polypropylene, these sturdy yet lightweight containers are 95kPa certified for air transport.

Leak-Resistant Dual Click-Tite™ Cap Seal

All containers feature fast screw threads with low-friction sealing surfaces for quick and easy container open and closure. Dual Click-Tite caps "click" twice when completely closed, indicating the container has been closed securely and is now leak-resistant. Inadvertent back-off of the cap during handling and transport is prevented by an anti-back-off system designed into the containers.

Freight Efficient Packaging

The 60mL x 53mm Dual Click-Tite containers are packed 560 per case in the same carton that holds 400 of the 90mL x 53mm containers. Choosing a 60mL x 53mm for your application not only reduces freight cost per container in this comparison – reduced container weight resulting from collecting and transporting less sample may also prove beneficial for inbound freight savings.

Fast, Accurate Temperature Reading

Containers are of a robust design while utilizing the absolute minimum of materials and energy during their production. Thinner sidewalls permit rapid temperature transfer for applications requiring immediate and accurate temperature sensing of the specimen.

Sterile and Tamper Evident

Sterile containers are complete with tamper evident label. Containers are stackable, graduated and include top and bottom knurls for ease of gripping with gloves.

FOR MORE INFORMATION

(909) 267-2000/(770) 614-2931

info@therapak.com

www.therapak.com



All containers are 53mm in diameter.

Part No.	Item	Packaging	Inner Pkg
740001	60mL Sterile White Cap w/ Label	560/cs	4 x 140
740002	60mL Sterile White Cap w/ Label and Temp	560/cs	4 x 140
740003	90mL Sterile White Cap w/ Label	400/cs	4 x 100
740004	90mL Sterile White Cap w/ Label and Temp	400/cs	4 x 100
740005	90mL Sterile White Cap w/ Label and Temp	400/cs	8 x 50
740006	120mL Sterile White Cap w/ Label	300/cs	4 x 75