


## Duramark VIP Shippers

The Duramark™ Vacuum Insulated Panel (VIP) Shipper is ideal for the transport of frozen (<-20°C) shipments.

This sturdy shipper has been drop tested under dry ice conditions and temperature tested to ISTA 7D standards. The innovative insulation panels require less dry ice while offering long lasting temperature control. There is a range of options to meet your needs: the small VIP shipper can maintain temperatures below -20°C for 2 days using only 3 pounds of dry ice while the 10.4 L size maintains -20°C for 6 days with 13 pounds of dry ice. The compact packaging requires less storage space, decreases freight costs and generates less waste.



3.2L VIP Shipper product performance with 4.7 pounds of Dry Ice at a constant 30°C ambient temperature.



Part No.	Item	ID	OD
562436	VIP Shipper - 3.2 Liter Capacity	8.4 x 4.4 x 5.3"	11.3 x 7.8 x 7.5"
562437	VIP Shipper - 8.2 Liter Capacity	9.1 x 8.6 x 6.4"	12.3 x 11.5 x 10"
562438	VIP Shipper - 10.4 Liter Capacity	9.1 x 8.6 x 8.1"	12.3 x 11.5 x 11.7"
562908	VIP Shipper - 18.5 Liter Capacity	11.3 x 8.9 x 10.4"	14 x 12.9 x 14.3"

### FOR MORE INFORMATION

(909) 267-2000/(770) 614-2931

info@therapak.com

www.therapak.com